

What's New in Version 3.5.8.26 10150 - April 2020

UPDATE METHOD – IMPORTANT	1
PREVIOUS UPDATES	1
PRINT AND SCAN MULTIPLE SCAN FORMS AT THE SAME TIME	2
Printing multiple forms	2
Scanning multiple forms	3
DARKER BOXES ON SCAN FORMS	4
ADDED AN ADDITIONAL PROFILE RESPONSE FOR SOMOS MILITARY SURVEY	4
ADDED OPTION TO SEPARATE SPINE IMPLANTS IN EXPORT	4
ADDED NEW KOOS AND HOOS 12 ITEM SCORES	5
Conversion of main HOOS/KOOS to 12 item scores	6
ADDED INTAKE SURVEY AND EXPORT FOR NJRR (AUSTRALIA) AND CRISTAL STUDY	6
New preop/intake registry surveys for hip, knee and shoulder arthroplasty	6
New custom export for NJRR data	7
ADDED DATE ICON AND POPULATE ADMISSION DATE TO SURGERY DATE	7
ADDED AN ADDITIONAL PAIN QUESTION TO VAS SANE AND PAIN SURVEYS	8
NEW FIELDS ADDED TO SURGERY DETAILS FOR HIP AND KNEE ARTHROPLASTY	8
Knee	8
Hip	9
PATIENT SUMMARY REPORT HAS FILTER ADDED FOR DETERIORATING SCORES	9
COMMUNICATION TO SYNCH DATABASE WITH WEB	10

UPDATE METHOD - IMPORTANT

This update requires the use of an .exe file for all users. You will find the link to download the .exe update file on the document with the Update instructions – this is in the 'Software Details and Update' page on the client section of the website.

PREVIOUS UPDATES


If you haven't updated recently and want to know what's new in previous updates this is on the client section of our website. If you don't yet have a log into this you can create one by registering on this link.

http://socratesortho.com/admin-registration




PRINT AND SCAN MULTIPLE SCAN FORMS AT THE SAME TIME

It is now possible to print the scan forms in a group rather than one at a time, and also to scan multiple scores/forms from only one screen, rather than scanning from each individual score screen.




Printing multiple forms

If you click on the print icon on the bottom of the screen window below will appear. Just select the scores you want, click print and they will all be printed out.




Scanning multiple forms


Click on the scan icon and the window below will open up. Click the scores that you want to scan, fill in the evaluation date or the follow-up period and click scan.



You will see a message on the screen asking you to put the first form into the scanner. They will be scanned in the order in which you select them.



When those pages have been scanned it will give you a message asking you to put the next form/s in the scanner, it will continue to do this until all the scores/forms you have selected have been scanned. If there are error messages from any of the forms you will get the usual messages, after it has finished scanning all the forms you will need to go into the screen of the score where there was an error and resolve whatever was causing it.




DARKER BOXES ON SCAN FORMS

Some of the scanned forms had lighter grey boxes than others and it was difficult for patients to see the boxes. The check boxes are now darker than they were before, you will need to reprint these from the new version.

ADDED AN ADDITIONAL PROFILE RESPONSE FOR SOMOS MILITARY SURVEY

Added an additional response for "not on profile"



ADDED OPTION TO SEPARATE SPINE IMPLANTS IN EXPORT


Spinal implants will now export in separate columns.

ADDED NEW KOOS AND HOOS 12 ITEM SCORES

Two new scores have been recently published which are now included in Socrates Hip and Knee Arthroplasty modules. The authors took 12 questions from the main versions of HOOS and KOOS. They are separated into 3 sub scores – Pain, Physical function and QOL with a total score out of 100.



Note that this is a NEW short score, it is not the KOOS/HOOS PS or JR both of which are already in Socrates.




The authors have concluded that "KOOS-12 and HOOS 12 was a reliable and valid alternative to KOOS/HOOS in TKR patients with moderate to severe OA"

The references and conclusions are on a separate attachment. No license is needed to use these scores.

Conversion of main HOOS/KOOS to 12 item scores.


If you would like to start using the new 12 item surveys in place of the full versions to reduce the burden on the patients it is possible convert your main HOOS/KOOS data and scores over to the new 12 item score. You **can only do this once,** it should be done when you have made a decision that going forward you will ONLY continue to use the 12 item scores. All the existing full HOOS and KOOS data will remain in the original screens but the 12 questions from those scores will be migrated at the same follow up time points to the new 12 item score where they will be rescored.

If you would like to do this, contact us for the password and instructions.



ADDED INTAKE SURVEY AND EXPORT FOR NJRR (AUSTRALIA) AND CRISTAL STUDY

New preop/intake registry surveys for hip, knee and shoulder arthroplasty A new survey as below has been added to the arthroplasty modules.






New custom export for NJRR data

This exports the PROMS completed by sites participating in the Australian joint registry pilot project. After exporting the surveys can be sent to the NJRR for importing to the registry. The PROMS are EQ5D 5L, OXFORD, KOOS 12 and HOOS 12, Cristal and the patient intake surveys.



ADDED DATE ICON AND POPULATE ADMISSION DATE TO SURGERY DATE

It's common for patients to be admitted the same day as surgery. For those sites that are using the admission date it is now possible to indicate this is the same day of surgery and the surgery date will also be automatically populated. The date icon to populate the current date has also been added to this field, it is in all modules.



ADDED AN ADDITIONAL PAIN QUESTION TO VAS SANE AND PAIN SURVEYS


Both the standard visual analogue pain score and the SANE/NORMAL surveys have had an additional question asking about pain at rest.




NEW FIELDS ADDED TO SURGERY DETAILS FOR HIP AND KNEE ARTHROPLASTY

Knee

Added new block site to this list




Added an additional field to record the type of incision lateral, central, medial, previous scar.




Hip

An additional field for direct superior approach has been added to the approach list.



Additional fields for releases added to the other release list. Anterior capsulotomy has been broken down into partial and complete, previous data entered under Anterior capsulotomy has been transferred to the new field - Complete anterior capsulotomy.



PATIENT SUMMARY REPORT HAS FILTER ADDED FOR DETERIORATING SCORES

The patient summary report now has a checkbox which enables you to generate the report only for patients who scores have deteriorated at any time point. If you are following your patients remotely this is a handy report to run from time to time, you can then recall any patients where you see a decline in their scores.



COMMUNICATION TO SYNCH DATABASE WITH WEB

The server currently gives a message when using a manual send to let users know that it is also communicating with the server to update any changes that have happened in the local database. For some users this process was slowing their system down. We have now changed this so it does this weekly and behind-the-scenes so you don't see the message.